

- Overview
- Monitoring Points
- Alerts 31
- Reports
- Calibration
- Audit Trail
- Manage
- Support

Anschutz Medical Campus / Circuit Management / CM\_AK32-143\_ULT:0401



CM\_AK32-143\_ULT:0401

Data Logger

19220700

How to Create a Notification on the DicksonOne Monitoring System

This will work for all notification -80 freezers, -20 freezers, Battery is low

Overview Alarms Settings

## Alarm Coverage

|                  | Sun                                   | Mon                                   | Tue                                   | Wed                                   | Thur                                  | Fri                                   | Sat                                   |
|------------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|
| Temperature CH:1 | is below -90.0°C<br>12:00AM – 11:59PM | is below -90.0°C<br>12:00AM – 11:59PM | is below -90.0°C<br>12:00AM – 11:59PM | is below -90.0°C<br>12:00AM – 11:59PM | is below -90.0°C<br>12:00AM – 11:59PM | is below -90.0°C<br>12:00AM – 11:59PM | is below -90.0°C<br>12:00AM – 11:59PM |
|                  | is above -65.0°C<br>12:00AM – 11:59PM | is above -65.0°C<br>12:00AM – 11:59PM | is above -65.0°C<br>12:00AM – 11:59PM | is above -65.0°C<br>12:00AM – 11:59PM | is above -65.0°C<br>12:00AM – 11:59PM | is above -65.0°C<br>12:00AM – 11:59PM | is above -65.0°C<br>12:00AM – 11:59PM |
| Battery CH:9     | is low<br>12:00AM – 11:59PM           | is low<br>12:00AM – 11:59PM           | is low<br>12:00AM – 11:59PM           | is low<br>12:00AM – 11:59PM           | is low<br>12:00AM – 11:59PM           | is low<br>12:00AM – 11:59PM           | is low<br>12:00AM – 11:59PM           |
| Data Logger      | not reporting<br>12:00AM – 11:59PM    | not reporting<br>12:00AM – 11:59PM    | not reporting<br>12:00AM – 11:59PM    | not reporting<br>12:00AM – 11:59PM    | not reporting<br>12:00AM – 11:59PM    | not reporting<br>12:00AM – 11:59PM    | not reporting<br>12:00AM – 11:59PM    |

## Applied Alarms (4)

+ Apply Alarm

+ Add a New Alarm

☐ Circuit Management\_logger not reporting

Edit

Pause

Remove

## When

Data Logger is not reporting

## Wait

5 readings

## Trigger

Warning

## Contact

Facilities Manager via Email  
<jeffrey.gulka@cuanschutz.edu>  
<Jason.radaelli@cuanschutz.edu>  
<heath.dejager@cuanschutz.edu>

## Schedule

24/7

## Channels (1)

Circuit Management - CM\_AK32-143\_ULT:0401 - Data Logger

Overview

Monitoring Points

Locations

Floor Plans

Equipment

Data Loggers

Channels

Alerts32

Reports

Calibration

Audit Trail

Manage

Support

Overview

Monitoring Points

Locations

Floor Plans

Equipment

Data Loggers

Channels

Alerts32

Reports

Calibration

Audit Trail

Manage

My Profile

Account

Users & Contacts

Locations

Data Loggers

Sensors

Gateways

Unassigned Data Loggers

Alarm Templates

Notification Policies

Integrations

Exports

Upload

API Keys

Billing

DICKSON

Environmental Monitoring + Compliance Experts

Privacy Terms of Use v1.12.2

Alarm Templates

Click on Add New Template

+ Add New Template

To get started, look for the "Overview" section on the left side of the screen. Click on "Manage." Once the screen opens with more options, select "Alarm Template." A new screen will appear; click on "Add New Template."

- 90 to -70 °C (9)
- Idrovo\_-85 to -75 C Pre Alarm (3)
- LEAD Center Not Reporting / Low Battery (17)
- LEAD Center -80C Alarm -90Cmin/-70Cmax (8)
- LEAD Center -20 C Alarm -35Cmin/-15Cmax (1)
- CTRC -65 - -90 (9)
- CTRC -25 - -15 (11)
- CTRC 2-8 (4)
- CTRC -15 - -25 prealarm (10)
- CTRC 2-8 pre-alarm (4)
- CTRC -65 - -90 pre alarm (10)
- Catenacci\_-70 to -90 C temp for Newest-80 freezer (1)
- Weinberg\_Battery alarm (18)
- Weinberg\_-150 Alarm (4)

Click on Manage

Click on Alarm Template

## Add Alarm Template

Name

Conditions

Notification Policy

Alarm Hours

Apply Channels

### 1) Name Your Template

Name \*

Test alarm for -80C Ultracold

Continue

When creating and naming your template, consider the following:

- Who will be using the template
- What assets or monitoring points these alarms are protecting

E.g., Vaccine Refrigerator High/Low Alarms

Create A Unique Name For Template

Click Continue

When creating and naming your template, consider the following:

Who will be using the template

What assets or monitoring points these alarms are protecting

E.g., Vaccine Refrigerator High/Low Alarms



## Add Alarm Template

Name

Conditions

Notification Policy

Alarm Hours

Apply Channels

### 2) Set Your Conditions

When

Select a Channel

Temperature

Relative Humidity

Differential Pressure

Dew Point

Battery

Data Logger

+ Add Another Condition

Continue

Select Temperature

Click on Continue

Configure the thresholds that should trigger alerts. You may have more than one condition on an alarm template. Specific notifications will be configured in the next step.

E.g., When temperature is above 0°C

Overview

Monitoring Points

Alerts 36

Reports

Calibration

Audit Trail

Manage

Support

## Add Alarm Template

Name

Conditions

Notification Policy

Alarm Hours

Apply Channels

## 2) Set Your Conditions

When

Temperature

is above

-65

°C

+ Add Another Condition

Continue

Back

Configure the thresholds that should trigger alerts. You may have more than one condition on an alarm template. Specific notifications will be configured in the next step.

E.g., When temperature is above 0°C

Make sure to set the temperature to Celsius.

The temperature should not be set Warmer than -65°C for a -80°C freezer.

There are two choices to be made

Is above or Is below

Click Add Another Condition

1. Choose Is above or Is below

2. Type in the temperature (don't forget the minus sign)


3. Change the setting to Celsius

4. Add Another Condition



- Overview
- Monitoring Points
- Alerts 34
- Reports
- Calibration
- Audit Trail
- Manage
- Support

# Add Alarm Template



## 2) Set Your Conditions

When

Temperature

is above

-65

°C

When

Temperature

is below

-85

°C

+ Add Another Condition

Continue Back

1. Choose Temperature

2. Is below

3. -85

4. Celsius

Click Continue


Configure the thresholds that should trigger alerts. You may have more than one condition on an alarm template. Specific notifications will be configured in the next step.

E.g., When temperature is above 0°C



- Overview
- Monitoring Points
- Alerts 34
- Reports
- Calibration
- Audit Trail
- Manage
- Support

## Add Alarm Template



### 3) Set Your Notification Policy

Use Existing

Create New

☐ Audible Alarm

+ Add Custom Escalation Group

Continue

Back

- Click on Create New
- Click on Add Custom Escalation Group
- Click Continue

A notification policy can have up to three escalations. Each escalation consists of three components: when to send the alert, the method to use, and who to contact. If you already have a policy configured, you can select that here, or you can create a new one.

E.g., After 5 readings, send an SMS to Contact.



# Add Alarm Template



## 3) Set Your Notification Policy

Use Existing   Create New

☐ Audible Alarm

### Escalation No. 01

Escalate as warning or excursion?

☒ Warning  
☐ Excursion

When should we let your contacts know?

Wait  readings before contacting them.

Who do you want to contact?

Notification No. 01   [Remove](#)

Select notification ▼

[+ Add Another Notification](#)

A notification policy can have up to three escalations. Each escalation consists of three components: when to send the alert, the method to use, and who to contact. If you already have a policy configured, you can select that here, or you can create a new one

E.g., After 5 readings, send an SMS to Contact.

Begin by creating a warning notification.

A reasonable wait time for the warning is four readings.

Then, select a type of notification.



# Add Alarm Template



## 3) Set Your Notification Policy

Use Existing   Create New

☐ Audible Alarm

### Escalation No. 01

Escalate as warning or excursion?

- ☒ Warning  
☐ Excursion

When should we let your contacts know?

Wait  readings before contacting them.

Who do you want to contact?

Select notification ^

Send email to

**Send SMS to**

Phone call to

Select contact or role ^

Roles

Biorepository Team

Primary Investigator

**Lab Manager**

Facilities Manager

Maintenance

Risk Management

Select contact method ^

**cell**

landline

work

phone

mobile

primary

home

This message is for setting up a text.

- 1. Select a notification – Send an SMS.
- 2. Choose a contact or role – select the one listed in your user account.
- 3. Choose a contact method – select cell phone.

Click continue

3) Set Your Notification Policy

Use Existing   Create New

☐ Audible Alarm

Escalation No. 01

Escalate as warning or excursion?

☒ Warning  
☐ Excursion

When should we let your contacts know?

Wait  readings before contacting them.

Who do you want to contact?

Notification No. 01   [Remove](#)

[+ Add Another Notification](#)

Escalation No. 02

Escalate as warning or excursion?

☐ Warning  
☒ Excursion

When should we let your contacts know?

Wait  readings before contacting them.

Who do you want to contact?

Notification No. 01   [Remove](#)

Send email to  
**Send SMS to**  
Phone call to

Roles  
Biorepository Team  
Primary Investigator  
**Lab Manager**  
Facilities Manager  
Maintenance

**cell**  
landline  
work  
phone  
mobile  
primary

[+ Add Custom Escalation Group](#)

[Continue](#)

[Back](#)

A notification policy can have up to three escalations. Each escalation consists of three components: when to send the alert, the method to use, and who to contact. If you already have a policy configured, you can select that here or create a new one.

E.g., After 5 readings, send an SMS to Contact.

1. Select a notification – Send an SMS.
2. Choose a contact or role – select the one listed in your user account.
3. Choose a contact method – select cell phone.

## Add Alarm Template

Name

Conditions

Notification Policy

Alarm Hours

Apply Channels

### 4) Set Alarm Schedule

Choose whether your alarms are scheduled 24/7 or only for specific days and times.

☒ Have this alarm active 24/7

Continue

[Back](#)

Here is where you can schedule when alarms will be received.

Click on Have this alarm active

Click Continue



- Overview
- Monitoring Points
- Alerts 29
- Reports
- Calibration
- Audit Trail
- Manage
- Support

# Add Alarm Template



## 4) Set Alarm Schedule

Choose whether your alarms are scheduled 24/7 or only for specific days and times.

☒ Have this alarm active 24/7

Select the day(s) you'd like this alarm active

Alarm schedules do not support transitions from one day to the next day. A start time and end time must be selected for each day that is added. If more shifts are needed, please click '+Add days and/or Times.'

☐ Sun ☐ Mon ☐ Tue ☐ Wed ☐ Thur ☐ Fri ☐ Sat

Active From:

12

:

00

AM

Until:

11

:

59

PM

+ Add Days and/or Times

☐ Enable Buffer

Applying a buffer to an alarm will allow alerts that begin outside of the scheduled window to trigger notifications immediately once the scheduled period begins. This setting can be used to escalate alerts between alarms and ensure that multiple shifts are covered.

☐ Enable After-Hours Notifications

Enabling after-hours notifications will allow back-to-normal alert notifications to be delivered to users when an alarm is triggered during the scheduled window, but is resolved outside of the alarm schedule. When after-hours notifications are not enabled, the system will not send notifications for alarms resolved outside of the alarm schedule.

Continue

Back

DICKSON  
Environmental Monitoring • Compliance Experts

Pick date and times as needed for schedule

- Overview
- Monitoring Points
- Alerts 29
- Reports
- Calibration
- Audit Trail
- Manage
- Support

# Add Alarm Template



## 5) Add Channels/Data Loggers

If you'd like to apply this template to the channels on a data loggers now, you can select from the list below or search to narrow down your selection. You can come back and do this at any time or apply a template to channels from a data logger's alarm page.

Search for data loggers, channels, or additional tags (ex. Fridge)

Q

Results (1611)

Applied to Template (0)

- ☐ P12-6th Floor - CORE LAB - Mercury Freezer 4 (-20) - 19122604 - Temperature CH:1

☐ P15\_11th Floor\_Weinberg Lab - (Bay 7 -20) Weinberg\_P15-11470C\_FRZ:0037 - 19210167 - Temperature CH:1

☐ P15\_11th Floor\_Weinberg Lab - (BGA -80) Weinberg\_P15-11470F\_ULT:0166 - 19209633 - Temperature CH:1

☐ P15\_11th Floor\_Weinberg Lab - (Butters -20) Weinberg\_P15-11380A\_FRZ:0023 - 19209631 - Temperature CH:1

☐ P15\_11th Floor\_Weinberg Lab - (Cartman -80) Weinberg\_P15-11390\_ULT:0078 - 19210427 - Temperature CH:1

☐ P15\_11th Floor\_Weinberg Lab - (Chef -20) Weinberg\_P15-11380B\_FRZ:0020 - 19209627 - Temperature CH:4

☐ L18\_10th Floor\_Johnson Lab - (Hall -20) L18-10th\_Johnson\_FRZ:0303 - 19122263 - Temperature CH:5

☐ P15\_11th Floor Weinberg Lab - (Henrietta -80) Weinberg\_P15-11370A\_ULT:0076 - 19210449 - Temperature

Selected results (left) will appear here and be applied to this template.

- Type in freezer names to find specific freezers
- OR
- Click on the box to assign the template to freezers