

Logger Not Reporting Data



⚠ No data within 24hr period.




****Logger Not Reporting Data****

1. Check the signal: No signal detected.
2. Turn off the logger.
3. Turn the logger back on and check for a WIFI signal.

Additionally, verify the sensor connection to ensure the connector is properly seated in the logger device.



WIFI Signal



Power Button



Sensor Connection



****Have you checked if your battery is low?****

When the battery runs out of power, all data stored in the logger will be lost.

It is the lab's responsibility to replace the battery as needed.

The logger requires two AA batteries for operation. Lab managers should set up a notification system for low battery levels to prevent any data loss.



Battery Power Level

Two AA Batteries

